

## Plan Element Nomenclature

---

Images .....	2
Beam Name (RS) / Field ID (MQ) .....	4
Beam Description (RS) / Field Name (MQ) .....	5
Beam Set .....	7
Plan .....	8
Points .....	9

Site Group: All	Original Date: December 2016
Issued By: RayStation Deployment Group	Revisions Date(s): None
Reviewed By: xxxx	Review Date:
Approved By: xxxx	Page 1

## Images

### Style

<Modality>\_<Technique>\_<Anatomic state/site>\_<Date>\_<External>

### *Modality*

Required. Indicates imaging modality.

Common modalities:

CT, MR, CBCT, PET

Examples:

CT\_07Sep2016, MR\_07Sep2016, CBCT\_07Sep2016, PET\_07Sep2016

### *Technique*

Optional. Indicates specific imaging technique.

Common techniques:

AVG, MIP, 0%, 50%, Arterial(?), Venous(?), Delayed(?), T1, T2, ADC, DCE, FDG

Examples:

MR\_T1\_07Sep2016, MR\_T2\_07Sep2016, MR\_ADC\_07Sep2016,

### *Anatomic State/Site*

Optional. Indicates the state or the site.

Common states/sites:

IBH, EBH, Preshift, Postshift, BladderFull, BladderEmpty, ERC

Examples:

CT\_IBH\_07Sep2016, CT\_EBH\_07Sep2016, CT\_Brain\_07Sep2016, CT\_Femur\_07Sep2016,  
CBCT\_Preshift\_07Sep2016, CBCT\_Postshift\_07Sep2016

Site Group: All	Original Date: December 2016
Issued By: RayStation Deployment Group	Revisions Date(s): None
Reviewed By: xxxx	Review Date:
Approved By: xxxx	Page 2

### ***Date***

Required. Indicates the date the image was acquired.

Examples:

CT\_06Aug2015, CT\_07Sep2016, CT\_08Sep2016

### ***External***

Required when present. Indicates if the image was acquired outside the Radiation Medicine Program.

CT\_07Sep2016\_External

### **Special Cases**

#### ***4DCT Groups***

Label as 4DCT\_<date>

Site Group: All	Original Date: December 2016
Issued By: RayStation Deployment Group	Revisions Date(s): None
Reviewed By: xxxx	Review Date:
Approved By: xxxx	Page 3

## Beam Name (RS) / Field ID (MQ)

### Style

<Beam Set Index><Beam Index><.Alternating Day Index>

### ***Beam Set Index***

Required. Indicates Beam Set and is consistent for a group of beams.

Examples:

A1, B1, C1

### ***Beam Index***

Required. Beam number within a Beam Set.

Examples:

A1, A2, A3

### ***Alternating Day Index***

Optional. Beam number within a Beam Set.

Examples:

A1.1, A2.1

Site Group: All	Original Date: December 2016
Issued By: RayStation Deployment Group	Revisions Date(s): None
Reviewed By: xxxx	Review Date:
Approved By: xxxx	Page 4

## Beam Description (RS) / Field Name (MQ)

### Style

<Beam Set Index>\_<Ref><PI><Bolus><G" Gantry Angle"> <Arc" Arc Span"" Arc Direction"> <C" Couch Angle"><Carriage Group>

### **Beam Set Index**

Required. Indicates Beam Set and is consistent for a group of beams.

Examples:

A\_G0, B\_Arc360CCW

### **Ref**

Optional. Indicates a beam is for reference only and not for treatment.

Examples:

A\_RefG0, B\_RefArc360CCW

### **Bolus**

Optional. Indicates a beam is to localize bolus only and not for treatment.

Examples:

A\_BolusG0, B\_BolusG90

### **PI**

Optional. Indicates a beam is for imaging only and not for treatment.

Examples:

A\_PIG0, B\_PIG270

Site Group: All	Original Date: December 2016
Issued By: RayStation Deployment Group	Revisions Date(s): None
Reviewed By: xxxx	Review Date:
Approved By: xxxx	Page 5

### ***Gantry Angle***

For 3D-CRT and IMRT fields. Indicates the gantry angle of the beam.

Examples:

A\_G0, B\_G320

### ***Couch Angle***

Optional. Not required when angle is zero. Indicates the couch angle of the beam.

Examples:

A\_G0C90, B\_G320C30

### ***Arc Span and Direction***

For VMAT fields. Indicates the length and direction of the arc.

Examples:

A\_Arc360CCW, B\_Arc90CW

### ***Multiple Treatment Fields at Same Gantry Angle***

Optional. For IMRT fields when carriage group splitting is required. Indicates apertures tend toward the Varian X1 jaw or X2 jaw.

Alternatively, can be used to differentiate multiple treatment fields that are geometrically similar but differ in energy.

Examples:

A\_G0X1, A\_G0X2

Site Group: All	Original Date: December 2016
Issued By: RayStation Deployment Group	Revisions Date(s): None
Reviewed By: xxxx	Review Date:
Approved By: xxxx	Page 6

## Beam Set

### Style

<Beam Set Index>.<Site><\_Laterality

### ***Beam Set Index***

Required. Indicates Beam Set and is consistent with the beam set index for the beams in the beam set.

Examples:

A.OPX\_L, B.Lung, C.TSpine

### ***Site***

Required. Indicates treated site, anatomic location, or volume. Site plus laterality must not exceed 14 characters.

Examples:

A.Breast\_IM\_L, B.Breast\_SA\_L

### ***Laterality***

Optional. Indicates left or right-sided treatment. Site plus laterality must not exceed 14 characters.

Examples:

A.Breast\_L, B.Breast\_R

Site Group: All	Original Date: December 2016
Issued By: RayStation Deployment Group	Revisions Date(s): None
Reviewed By: xxxx	Review Date:
Approved By: xxxx	Page 7

## Plan

### Style

<Treatment>\_<Beam Set Indices>

### ***Treatment***

Required. Indicates the site and type of treatment.

Examples:

Breast\_L\_A, Breast\_R\_B

### ***Beam Set Indices***

Required. Indicates the Beam Sets contained within the plan and is consistent with the beam set index for the beams in the beam sets.

Examples:

Prostate \_A, ProstateNodes\_AB

Site Group: All	Original Date: December 2016
Issued By: RayStation Deployment Group	Revisions Date(s): None
Reviewed By: xxxx	Review Date:
Approved By: xxxx	Page 8



## Points

### Laser Position

The point that represents the intersection of the lasers at CT-simulation is always named CT SETUP

### Prescription

Consists of “ICRU” and a letter indicating the associated beams

#### *Style*

ICRU\_<Beam Set Index>

#### *Beam Set Index*

Required. Indicates Beam Set and is consistent with the beam set index for the beams in the beam set.

Examples:

ICRU\_A, ICRU\_B

### Isocentres

Consists of “Iso”, a letter indicating the associated beams, and an optional description

#### *Style*

Iso\_<Beam Set Index>\_<description>

#### *Beam Set Index*

Required. Indicates Beam Set and is consistent with the beam set index for the beams in the beam set.

Examples:

Iso\_A, Iso\_B

#### *Description*

Optional. Used when an additional description is required to distinguish the isocentres.

Examples:

Iso\_A\_Sup, Iso\_B\_Inf

Site Group: All	Original Date: December 2016
Issued By: RayStation Deployment Group	Revisions Date(s): None
Reviewed By: xxxx	Review Date:
Approved By: xxxx	Page 9

### Dose Specification Points

Consists of “DSP” and a description

#### *Style*

**DSP\_<description>**

#### *Description*

Required. The description can refer to the location of a prescription point, an alternative point, or a specific associated beam.

Examples:

DSP\_ICRU\_A, DSP\_Iso\_A, DSP\_G120

### Others (markers, tattoos, automated points)

Points for tattoos and markers refer to location on body that they indicate.

“Ref” will prefix points created for reference or measurements.

“Reg” will prefix points created for image registration.

Examples

INF, SUP, LAT, MED, ANT

RefTTH, RefCouch

RegAnterior, Reg1

Site Group: All	Original Date: December 2016
Issued By: RayStation Deployment Group	Revisions Date(s): None
Reviewed By: xxxx	Review Date:
Approved By: xxxx	Page 10